



info@eqas.gr

## COLOMBO - Automatic Coagulation Analyzer

COLOMBO is a **fully automated coagulation analyzer**, with photo-optical clot detection.

COLOMBO has high sensitive photometer designed to perform all routine test such as **PT, APTT, TT, Fibrinogen**, with a light source provided by a solid state three wave-length emitter at 460 nm, 510 nm, 630 nm which ensure accurate results.

- 2 Measuring Channel: 1 Clotting, 1 Chromogenic
- Parameters: PT, APTT, TT, FIB
- Photo-optical clot detection system for a real walk away test performing
- High performances graphic software running under windows xp
- Quality controls,databases,graphical results.

