

# CATALOGUE

**ScreenItalia**  
Test per la salute

# SUMMARY

DRUG TESTS 2	DISPOSABLE ALCOHOL TESTS 29	ALCOHOL TESTERS 31	FLUORESCENCE TESTS 33	PREGNANCY/ GYNECOLOGY 43	INFECTIVE DISEASES 46
66 INFLAMMATORY MARKERS	68 TUMORAL MARKERS	71 CARDIAC MARKERS	73 URINE ANALYSIS	76 OTHER TESTS	78 RAPID TEST READERS



...DRUG...

...TESTS...

2

DRUG TESTS

# SCREEN DROGA TEST CAPELLO

(MULTI DRUG TEST HAIR)

Product Code: SC-1170-20

A rapid test for the qualitative detection of Cotinine, d-Methamphetamine, Morphine, Cocaine in human hair. Hair detection has unique advantages in the field of drug testing. For professional in vitro diagnostic use only.



# SCREEN DROGA TEST SUPERFICI/POLVERI

(MULTI DRUG RAPID TEST SURFACE SOLIDS)

Product code: SC-0371-25


Screen Surface/Solids Drug Test is a rapid test for surfaces or solids that can be performed without the use of an instrument. The test utilizes monoclonal antibodies to selectively detect elevated levels of specific drugs on surfaces and solids. Test is intended to be used as analytical device to detect drugs on surfaces.



# SCREEN TEST NICOTINA

(COTININE RAPID TEST PANEL)

Product Code: SC-1019-40



A rapid test for the qualitative detection of Cotinine (Nicotine metabolite). For determination of smoking status only. The COT Rapid Test Cassette is a rapid chromatographic immunoassay for the detection of Cotinine. Not intended for medical diagnostic use.

# SCREEN TEST ADULTERAZIONE

(ADULTERATION STRIPS)

Product code: SC-0807-25




The test allows verifying whether the urine sample has been counterfeited through Creatinine, pH, Nitrites, Glutaraldehyde, specific gravity (viscosity of the urine sample), Oxidants and Bleach. Each plastic strip contains n. 6 chemical reacting swabs. For professional in vitro diagnostic use only.

# SCREEN 1 BENZODIAZEPINE

(ONE STEP DRUG TEST CARD)

Product Code: SC-0036-25

Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.



# SCREEN 1 BUPRENORFINE

(ONE STEP DRUG TEST CARD)

Product Code: SC-0029-25

Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.





# SCREEN 1 KETAMINA

(ONE STEP DRUG TEST CARD)

Product code: SC-0067-25



Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.

# SCREEN 1 MARIJUANA

(ONE STEP DRUG TEST CARD)

Product code: SC-0098-25



Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.

# SCREEN 1 AMFETAMINA

(ONE STEP DRUG TEST CARD)

Product code: SC-1781-25

Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.



# SCREEN 1 BARBITURICI

(ONE STEP DRUG TEST CARD)

Product code: SC-1774-25

Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.



# SCREEN 1 METAMFETAMINA

(ONE STEP DRUG TEST CARD)

Product code: SC-1767-25



Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.

# SCREEN 1 COCAINA

(ONE STEP DRUG TEST CARD )

Product code: SC-0043-25



Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.




# SCREEN 1 MDMA/EXSTASY

(ONE STEP DRUG TEST CARD)

Product Code: SC-0487-25

Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.




# SCREEN 1 METADONE

(ONE STEP DRUG TEST CARD)

Product Code: SC-0074-25

Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.





# SCREEN 1 OPIACEI

(ONE STEP DRUG TEST CARD)

Product Code: SC-0081-25



Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.

# SCREEN 1 K2-SPICE

(ONE STEP DRUG TEST CARD)

Product code: SC-0050-25




Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.

# SCREEN 1 FENTANIL

(ONE STEP DRUG TEST CARD)

Product Code: SC-1828-25

Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.




# SCREEN 1 PROPOSSIFENE

(ONE STEP DRUG TEST CARD)

Product Code: SC-2030-25

Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.



# SCREEN 1 TRAMADOLO

(ONE STEP DRUG TEST CARD)

Product Code: SC-2047-25



Rapid test for the qualitative detection of drug in human urine. This single test is available to detect: Cocaine and metabolites, Buprenorphine, Amphetamine, Methamphetamine (including Ecstasy MDMA), Ketamine (KET), Marijuana (including Cannabinoids THC), Benzodiazepines (BZO), Oxycodone (OXY), Opiates (including Heroin, Morphine and metabolites), etc. This test can be customized based on the substance required.

# SCREEN 5

(MULTI LINE DRUG RAPID TEST CASSETTE)

Product code: SC-0128-25



The test detects 5 substances simultaneously, result is specific to each substance analyzed. Detects: Propoxyphene, Ketamine, Fentanyl, Buprenorphine, Tricyclic Antidepressants. For professional in vitro diagnostic use only.



# SCREEN 5

(MULTI DRUG RAPID TEST PANEL)

Product code: SC-0104-25

The test detects 5 substances simultaneously, result is specific to each substance analyzed. Detects: Cocaine and metabolites, Methamphetamine (included Ecstasy MDMA), Amphetamine, Marijuana, Opiates. For professional in vitro diagnostic use only.



# SCREEN 5

(MULTI LINE DRUG RAPID TEST CASSETTE )

Product code: SC-0111-25


The test detects 5 substances simultaneously, result is specific to each substance analyzed. Detects: Oxycodone, EDDP, Tramadol, Phencyclidine, Synthetic Marijuana K2. For professional in vitro diagnostic use only.



# SCREEN 7+AD

(MULTI DRUG RAPID TEST PANEL)

Product code: SC-0135-25 / SC-0821-25



The test detects 7 substances simultaneously, results for each substance. Detects: Cocaine (COC 300), Amphetamine (AMP 500), Methamphetamine (MET 500), Marijuana (THC 50), Opiates (MOP 300), Methylenedioxymethamphetamine (MDMA 500), Methadone (MTD 300) includes control for Creatinine e pH (adulteration). For professional in vitro diagnostic use only.

# SCREEN CUP 7+AD

(MULTI DRUG RAPID TEST CUP)

Product code: SC-0142-25




The test detects 7 substances simultaneously plus adulterations: Cocaine (COC 300), Amphetamine (AMP 500), Methamphetamine (MET 500), Marijuana (THC 50), Opiates (MOP 300), Methylenedioxymethamphetamine (MDMA 500), Methadone (MTD 300). Includes control for Creatinine e pH (adulteration). For professional in vitro diagnostic use only.

# SCREEN CUP 7+ADU

(MULTI DRUG RAPID TEST CUP)

Product code: SC-1224-25

Urine Drug Test Cup with Activation Key 7 + AD. The test detects 7 substances simultaneously plus 6 adulterants: Cocaine and metabolites (COC 300), Amphetamine (AMP 500), Methamphetamine (MET 500), Marijuana (Including THC Cannabinoids) (THC 50), Opiates (Including Heroin, Morphine and metabolites) (MOP 300), Methylenedioxymethamphetamine (Ecstasy, MDMA 500), Methadone (MTD 300), also checks the validity of the urine sample through an integrated adulteration test: OX, SG, pH, NIT, GLUT and CRE. For professional in vitro diagnostic use only.



# SCREEN 8

(MULTI LINE DRUG RAPID TEST CASSETTE)

Product code: SC-0173-25

The test detects 8 substances simultaneously plus adulterations: Cocaine (COC 300), Amphetamine (AMP 500), Methamphetamine (MET 500), Buprenorphine (BUP 10), Marijuana (THC 50), Opiates (MOP 300), Methylenedioxymethamphetamine (Ecstasy, MDMA 500), Methadone (MTD 300). For professional in vitro diagnostic use only.





# SCREEN 8

(MULTI LINE DRUG RAPID TEST CASSETTE)

Product code: SC-0166-25



The test detects 8 substances simultaneously plus adulterations: Cocaine (COC 300), Amphetamine (AMP 500), Methamphetamine (MET 500), Buprenorphine (BUP 10), Marijuana (THC 50), Opiates (MOP 300), Methylenedioxymethamphetamine (Ecstasy, MDMA 500), Methadone (MTD 300). For professional in vitro diagnostic use only.

# SCREEN 8

(MULTI DRUG RAPID TEST PANEL)

Product code: SC-0180-25



The test detects 8 substances simultaneously: Cocaine, Amphetamine, Marijuana, Opiates, Benzodiazepine, Barbiturates, Methadone and Tricyclic Antidepressants. For professional in vitro diagnostic use only.

# SCREEN 8+AD

(MULTI DRUG RAPID TEST PANEL)

Product code: SC-0418-25

The test detects 8 substances simultaneously plus Creatinine and pH: Cocaine and metabolites (COC 300), Amphetamine (AMP 500), Methamphetamine (MET 500), Buprenorphine and metabolites (BUP 5), Marijuana (including Cannabinoids THC) (THC 50), Opiates (including Heroin, Morphine and metabolites) (MOP 300), Methylenedioxymethamphetamine Ecstasy (MDMA 500), Methadone (MTD 300). For professional in vitro diagnostic use only.



# SCREEN 8+AD

(MULTI DRUG RAPID TEST PANEL)

Product code: SC-0159-25 / SC-0838-25

The test detects 8 substances simultaneously plus Creatinine and pH: Cocaine (COC 300), Amphetamine (AMP 500), Methamphetamine (MET 500), Buprenorphine (BUP 10), Marijuana (THC 50), Opiates (MOP 300), Methylenedioxymethamphetamine (Ecstasy, MDMA 500), Methadone (MTD 300). For professional in vitro diagnostic use only.





# SCREEN CUP 8+AD

(MULTI DRUG RAPID TEST CUP)

Product code: SC-0203-25 / SC-0814-25




The test detects 8 substances simultaneously plus adulterations: Cocaine (COC 300), Amphetamine (AMP 500), Methamphetamine (MET 500), Buprenorphine (BUP 10), Marijuana (THC 50), Opiates (MOP 300), Methylenedioxymethamphetamine (Ecstasy, MDMA 500), Methadone (MTD 300) and includes control for Adulterants. Adulteration test is based on the color response of chemical indicators in the presence of Creatinine (CR), pH. For professional in vitro diagnostic use only.

# SCREEN CUP 8+A+AD

(MULTI DRUG RAPID TEST CUP)

Product code: SC-0197-25



The test detects 8 substances simultaneously plus adulterations: Cocaine (COC 300), Amphetamine (AMP 500), Methamphetamine (MET 500), Buprenorphine (BUP 10), Marijuana (THC 50), Opiates (MOP 300), Methylenedioxymethamphetamine (Ecstasy, MDMA 500), Methadone (MTD 300), includes control for Alcohol. Adulteration test is based on the color response of chemical indicators in the presence of adulterants. Creatinine (CR), pH and Alcohol test ETG. For professional in vitro diagnostic use only.

# SCREEN 9

(MULTI LINE DRUG RAPID TEST CASSETTE )

Product code: SC-0210-25

The test detects 9 substances simultaneously plus adulterations: Methadone, Morphine, Cocaine, Marijuana (THC), Amphetamine, Benzodiazepine, Tricyclic Antidepressants, Barbiturates, Methylenedioxymethamphetamine (MDMA). For professional in vitro diagnostic use only.



# SCREEN 10

(MULTI LINE DRUG RAPID TEST CASSETTE)

Product code: SC-0258-25

The test detects 10 substances simultaneously: Marijuana (THC 50ng/ml), Cocaine and metabolites (COC 300ng/ml), Amphetamine (AMP 500ng/ml), Methamphetamine (mAMP 500ng/ml), Opiates/Morphine (300ng/ml), Buprenorphine (BUP 10ng/ml), Benzodiazepines (BZO 300ng/ml), Methadone (MTD 300ng/ml), Ecstasy (MDMA 500ng/ml), Tricyclic antidepressants (TCA 100ng/ml). For professional in vitro diagnostic use only.



# SCREEN 10

(MULTI DRUG RAPID TEST PANEL)

Product Code: SC-0234-25



The test detects 10 substances simultaneously: Cocaine and metabolites (300), Amphetamines (1000), Methamphetamines (Ecstasy MDMA) (1000), Marijuana (Cannabinoids THC) (50), Opiates (Morphine, Heroin and metabolites) (2000), Buprenorphine (10), Benzodiazepines (300), Barbiturates (300), Methadone (300), Tricyclic antidepressants (1000). For professional in vitro diagnostic use only.

# SCREEN 10

(MULTI LINE DRUG RAPID TEST CASSETTE)

Product code: SC-0241-25



The test detects 10 substances simultaneously: Marijuana (THC 50ng/ml), Cocaine and metabolites (COC 300ng/ml), Amphetamine (AMP 1000ng/ml), Methamphetamine (mAMP 500ng/ml), Opiates/Morphine (MOP300ng/ml), Phencyclidine (PCP 25ng/ml), Benzodiazepines (BZO 300ng/ml), Barbiturates (BAR 300ng/ml), Methadone (MTD 300ng/ml), Ecstasy (MDMA 500ng/ml). For professional in vitro diagnostic use only.



# SCREEN 10

(MULTI LINE DRUG RAPID TEST CASSETTE)

Product code: SC-0265-25

The test detects 10 substances simultaneously: Buprenorphine, Phencyclidine, Tricyclic Antidepressants, Tramadol, Ketamine, EDDP, Fentanyl, K2/Spice, Propoxyphene, Oxycodone. For professional in vitro diagnostic use only.




# SCREEN 10

(MULTI DRUG RAPID TEST PANEL)

Product Code: SC-0227-25

The test detects 10 substances simultaneously: Cocaine and metabolites, Amphetamine, Methamphetamine (including Ecstasy MDMA), Marijuana, Opiates, Benzodiazepine, Phencyclidine, Barbiturates, Methadone and Tricyclic antidepressants. For professional in vitro diagnostic use only.



# SCREEN CUP 10+A+AD

(MULTI DRUG RAPID TEST CUP)

Product code: SC-0272-25



The test detects 10 substances simultaneously plus adulterations and Alcohol test ETG: Amphetamine, Barbiturates, Buprenorphine, Benzodiazepine, Cocaine, MDMA (Ecstasy), Methadone, Opiates, Phencyclidine, Thc (Marijuana). For professional in vitro diagnostic use only.

# SCREEN CUP 10+A+AD

(MULTI DRUG RAPID TEST CUP )

Product code: SC-0289-25



The test detects 10 substances simultaneously plus adulterations and Alcohol test ETG: Amphetamine, Methamphetamine, Barbiturates, Buprenorphine, Benzodiazepine, Cocaine, Opiates, Methadone, Ketamine, Thc (Marijuana). For professional in vitro diagnostic use only.

# SCREEN 12

(MULTI LINE DRUG RAPID TEST CASSETTE )

Product code: SC-0319-25

The test detects 12 substances simultaneously: Marijuana (THC 50ng/ml), Cocaine and metabolites (COC 300ng/ml), Oxycodone (OXY 100ng/ml), Methamphetamine (mAMP 500ng/ml), Opiates/Morphine (MOP300ng/ml), Phencyclidine (PCP 25ng/ml), Benzodiazepines (BZO 300ng/ml), Methadone (MTD 300ng/ml), Ecstasy (MDMA 500ng/ml), Ketamine (KET 500ng/m), Tramadol (TML 100ng/ml), K2/Spice (K2 50ng/ml). For professional in vitro diagnostic use only.



# SCREEN 12

(MULTI DRUG RAPID TEST PANEL)

Product code: SC-0326-25

The test detects 12 substances simultaneously: Opiates/Morphine (MOP300ng/ml), Ecstasy (MDMA 500ng/ml), Oxycodone (OXY 100ng/ml), Phencyclidine (PCP 25ng/ml), Cocaine and metabolites (COC 300ng/ml), Methadone (MTD 300ng/ml), Barbiturates (BAR 300ng/ml), Benzodiazepines (BZO 300ng/ml), Amphetamine (AMP 1000ng/ml), Methamphetamine (MET 500ng/ml), Tricyclic antidepressants (TCA 1000ng/ml), Marijuana (THC 50ng/ml). For professional in vitro diagnostic use only.





# SCREEN CUP 12+A+AD

(MULTI DRUG RAPID TEST CUP)

Product code: SC-0296-25




The test detects 12 substances simultaneously plus adulterations and Alcohol test ETG: Barbiturates, Buprenorphine, Opiates, Benzodiazepine, Cocaine, Amphetamine, Methamphetamine, Methadone, Oxycodone, Ketamine, Propoxyphene, Thc (Marijuana). For professional in vitro diagnostic use only.

# SCREEN CUP 12+A+AD

(MULTI DRUG RAPID TEST CUP)

Product code: SC-0302-25



The test detects 12 substances simultaneously plus adulterations and Alcohol test ETG: Barbiturates, Buprenorphine, Cocaine, Benzodiazepines, MDMA (Ecstasy), Methamphetamine, Opiates, Methadone, Oxycodone, Phencyclidine, Propoxyphene, Thc (Marijuana). For professional in vitro diagnostic use only.

# SCREEN 6-5

(MULTI DRUG RAPID TEST ORAL FLUID)

Product Code: SC-0388-25

Saliva drug test. The test detects 6 substances: Opiates and metabolites, Cocaine and metabolites, Cannabinoids Marijuana (THC), Amphetamine, Methylenedioxymethamphetamine (MDMA) and Methadone. For professional in vitro diagnostic use only.



# SCREEN MIDSTREAM 6-5

(MULTI DRUG RAPID TEST ORAL FLUID)

Product code: SC-0449-25

A rapid test for the simultaneous, qualitative detection of multiple drugs and drug metabolites in human saliva. The Multi-Drug Rapid Test for Opiates (OPI40), Cocaine (COC30), Marijuana (THC25), Amphetamine (AMP40), Methamphetamine (MET40), Oxycodone (OXY40), Methylenedioxymethamphetamine (MDMA50). For professional in vitro diagnostic use only.





# SCREEN MIDSTREAM 6-5-MED

(MULTI DRUG RAPID TEST ORAL FLUID )

Product code: SC-0333-25



A rapid test for the simultaneous, qualitative detection of multiple drugs and drug metabolites in human saliva. The Multi-Drug Rapid Test for Opiates, Cocaine, THC/Marijuana, Amphetamine, Methamphetamine, MDMA. For professional in vitro diagnostic use only.

# SCREEN 6-52

(MULTI DRUG RAPID TEST ORAL FLUID)

Product Code: SC-0395-25



Saliva drug test. The test detects 6 substances: Opiates, Oxycodone, Cocaine, Cannabinoids Marijuana (THC), Methamphetamine including Methylenedioxymethamphetamine (MDMA, Ecstasy), Amphetamines. For professional in vitro diagnostic use only.

# SCREEN MIDSTREAM 6-52

(MULTI DRUG RAPID TEST ORAL FLUID)

Product code: SC-0364-25

A rapid test for the simultaneous, qualitative detection of multiple drugs and drug metabolites in human saliva. The Multi-Drug Rapid Test for Opiates, Cocaine, Thc/Marijuana, Amphetamine, Methamphetamine/MDMA, Synthetic Marijuana K2. For professional in vitro diagnostic use only.



# SCREEN MIDSTREAM 6-53

(MULTI DRUG RAPID TEST ORAL FLUID)

Product code: SC-0357-25

A rapid test for the simultaneous, qualitative detection of multiple drugs and drug metabolites in human saliva. The Multi-Drug Rapid Test for Opiates, Cocaine, Marijuana (THC), Amphetamine, Methamphetamine/MDMA, Benzodiazepines. For professional in vitro diagnostic use only.



# SCREEN MIDSTREAM 6-54

(MULTI DRUG RAPID TEST ORAL FLUID)

Product code: SC-0340-25



A rapid test for the simultaneous, qualitative detection of multiple drugs and drug metabolites in human saliva. The Multi-Drug Rapid Test for Opiates, Cocaine, Marijuana (THC), Amphetamine, Methamphetamine/MDMA, Methadone. For professional in vitro diagnostic use only.

# SCREEN CUP 10-S

(MULTI DRUG RAPID TEST ORAL FLUID)

Product Code: SC-0005-25



Screen Cup 10-S is a lateral flow chromatographic immunoassay for the qualitative detection of multiple drugs and drug metabolites in saliva. The Multi-Drug Rapid Test is a rapid, saliva screening test that can be performed without the use of an instrument. The test utilizes monoclonal antibodies to selectively detect elevated levels of specific drugs in human saliva. The Multi-Drug Rapid Test for Opiates, Methadone, Cocaine, Thc/Marijuana, Amphetamine, Metamphetamine, Mdma, K2 Synthetic Cannabinoid, Ketamine, Alcohol.

# ALCOHOL TESTS...

**29**

**DISPOSABLE ALCOHOL TESTS**



# SCREEN ALCOL TEST ETG

(ETHYL GLUCURONIDE TEST PANEL)

Product code: SC-0012-40



The Ethyl Glucuronide (ETG) Rapid Test Panel (Urine) is a rapid chromatographic immunoassay for the detection of Ethyl Glucuronide in human urine. Ethyl Glucuronide (ETG) is a metabolite of ethyl alcohol which is formed in the body by glucuronidation following exposure to ethanol, such as by drinking alcoholic beverages. ETG is a more accurate indicator of the recent exposure to alcohol than measuring for the presence of ethanol itself as it can be measured in urine by traditional laboratory methods (GC/MS or LC/MS) up to approximately 80 hours after ethanol is ingested (greatly dependent on the amount of alcoholic and non-alcoholic drinks). This assay provides only a preliminary analytical test result. For professional in vitro diagnostic use only.

# ETIL PROFESSIONAL

(SALIVA ALCOHOL TEST STRIP)

Product code: SC-1002-40



The Screen Etil Professional rapid test is an chemical method test for the rapid and sensitive detection of alcohol (alc) (ethyl alcohol/ethanol) in saliva and other liquid specimen. By examination of saliva specimen you can estimate the blood alcohol concentration. This chemical method gives a qualitative result, i.e. the test shows if specimens e.g. saliva contain alcohol or not. The assay detects ethanol at concentration of 0.1‰ or 10 mg/dL respectively by colour change of the reaction or test pad. The blue colour intensity increases with the alcohol concentration in the sample. Thus, Etil Pro CE produces a colour change in the presence of alcohol in specimen ranging from light yellow at negative specimens via a light blue, middle blue to a dark bluegrey colour at high alcohol concentration.

# ALCOHOL TESTERS

**31** | ALCOHOL TESTERS

# BAC 100 ALCOHOL DETECTOR

(ALCOHOL BREATH ANALYZER)

Product code: WTO5020102



This equipment analyzes the exhaled alveolar air, both in the active and passive mode (i.e., if the patient is unconscious). Product available for: Public Officers, Industrial Medicine, Hospitals, Analysis, Labs, Emergency Departments, etc. Professional electronic breath analyzer with incorporated printer. The breath analyzer is presently employed for professional use both by the industrial medicine and by police officers, as a precursor. Its precision and integrated thermal printer make it an extremely versatile and useful portable instrument. It comes with a case including: the breath analyzer, 6 rolls of thermal paper, 10 mouthpieces (individually wrapped), one passive mode mouthpiece, battery charger for automobile, electric current battery charger (for 8 rechargeable batteries).

# ALCOTRUE M

(ALCOHOL BREATH ANALYZER)

Product code: 5040112002 / 4040112007 Thermal Paper Printer



Portable Breathalyzer with electrochemical sensors of class I with a margin of error of  $\pm 2\%$ , printer and disposable mouthpiece. The breathalyzer has been realized with Medical CE according to CEE 93/42 on medical devices (MDD 93/42), specific for the Occupational Medicine is certified as a medical device in Class I with a measuring function. Join the Directory Medical Devices of the Ministry of Health in Italy with the Progressive System: 1126581, GMDN 36041 nomenclature, classification CND W0201069003. Its accuracy and the thermal print make the instrument extremely versatile and convenient to carry. It comes in a carrying case with inside: The Breathalyzer, 25 mouthpieces packaged individually and 4 rechargeable batteries.



# FLUORESCENCE ...TESTS...

33

FLUORESCENCE TESTS



# SCREEN IFA TEST COVID-19 IGG/IGM

(FLUORESCENCE IMMUNOASSAY COVID-19 IGG/IGM)

Product code: SC-1309-20



The Covid-19 IgG/IgM Test Cassette is a Fluorescence Immunoassay for the qualitative detection of IgG and IgM antibodies to Sars-CoV-2 in human whole blood, serum or plasma with the use of Fluorescence Immunoassay Analyzer as an aid in the diagnosis of Covid-19 infection. The Covid-19 IgG/IgM Test is an aid in the diagnosis of patients with suspected Sars-CoV-2 infection in conjunction with clinical presentation and the results of other laboratory tests. Results from the Covid-19 IgG/IgM Test should not be used as the sole basis for diagnosis. For professional in vitro diagnostic use only.

# SCREEN IFA TEST COVID-19 AG

(FLUORESCENCE IMMUNOASSAY COVID-19 ANTIGEN)

Product code: SC-1293-20



The Covid-19 Antigen Test Cassette is a Fluorescence Immunoassay for the qualitative detection of Sars-CoV-2 antigens present in human rhino pharynx with the use of Fluorescence Immunoassay Analyzer. The Covid-19 Antigen Test Cassette (Nasopharyngeal Swab) is based on fluorescence immunoassay for the qualitative detection of Sars-CoV-2 antigens in nasopharyngeal swab specimens from individuals with suspected Sars-CoV-2 infection in conjunction with clinical presentation and the results of other laboratory tests. For professional in vitro diagnostic use only.

# SCREEN IFA ANALYZER SMART

(FLUORESCENCE IMMUNOASSAY ANALYZER)

Product code: SC-1323-01

Fluorescence Immunoassay Analyzer is an analyzer that based on detection of fluorescence emitted during an immunoassay with antigen-antibody interaction. Fluorescence Immunoassay Analyzer offers the advantages of high accuracy, strong stability and fast results. Screen IFA Analyzer Fluorescence Immunoassay Analyzer should be used only with in vitro diagnostic tests manufactured by Screen Italia Srl as per package insert provided with specific test kits used.



# SCREEN IFA ANALYZER

(FLUORESCENCE IMMUNOASSAY ANALYZER)

Product code: SC-1491-01

Fluorescence Immunoassay Analyzer is an analyzer that based on detection of fluorescence emitted during an immunoassay with antigen-antibody interaction. It is designed to provide quantitative or qualitative test results by the examination of human samples with specific in vitro diagnostic test units including Inflammation markers, tumor markers, nephrology, diabetes, cardiac markers, coagulation, endocrinology, autoimmunity, infectious diseases etc. Offers the advantages of high accuracy, strong stability and fast results. Fluorescence Immunoassay Analyzer should be used only with in vitro diagnostic tests manufactured by Screen Italia Srl.





# SCREEN IFA TEST CK-MB

(FLUORESCENCE IMMUNOASSAY CK-MB)

Product code: SC-1477-20



A test for the diagnosis of Creatine Kinase MB (CK-MB) by measuring Creatine Kinase MB in whole blood, serum or plasma with the use of Fluorescence Immunoassay Analyzer. For professional in vitro diagnostic use only. The CK-MB Test Cassette (whole blood/serum/plasma) is a chromatographic immunoassay for the quantitative detection of human CK-MB in whole blood, serum or plasma as an aid in the diagnosis of Myocardial Infarction (MI). For professional in vitro diagnostic use only.

# SCREEN IFA TEST CRP

(FLUORESCENCE IMMUNOASSAY C-REACTIVE PROTEIN)

Product code: SC-1361-20



A Fluorescence Immunoassay test kit for the diagnosis of inflammatory condition and ACS by measuring CRP/hs-CRP in whole blood, serum or plasma with the use of Fluorescence Immunoassay Analyzer. The CRP Test Cassette (whole blood/serum/plasma) is based on Fluorescence Immunoassay for the quantitative determination of C-Reactive Protein (CRP) in serum, plasma or whole blood as an aid in the evaluation of infection, tissue injury and inflammatory disorders along with measurement of high sensitivity CRP (hs-CRP) for evaluation of acute coronary syndromes (ACS). For professional in vitro diagnostic use only.



# SCREEN IFA TEST MICROALBUMINURIA

(FLUORESCENCE IMMUNOASSAY MICROALBUMIN)

Product code: SC-1330-20

A Fluorescence Immunoassay test kit for measuring Microalbumin in urine with the use of Fluorescence Immunoassay Analyzer. The steady expulsion of small quantities of albumin with the urine can be the first sign of kidney damage. In the healthy kidney, Albumin is usually glomerular filtrated and tubular reabsorbed so that it is hardly detectable in urine. With a damaged kidney this process is disordered. The expulsion of Albumin in the range of 20– 200 mg/L is characterized as Microalbuminuria. For professional in vitro diagnostic use only.




# SCREEN IFA TEST HBA1C

(FLUORESCENCE IMMUNOASSAY HBA1C)

Product code: SC-1347-20

A Fluorescence Immunoassay test kit for measuring HbA1c value in whole blood with the use of Fluorescence Immunoassay Analyzer. The measure of HbA1c is recommended as a marker of long-term metabolic control in persons with diabetes mellitus. This test can be used as an aid in the diagnosis of diabetes and as an aid in identifying patients who may be at risk for developing diabetes. For professional in vitro diagnostic use only.



# SCREEN IFA TEST LH

(FLUORESCENCE IMMUNOASSAY HORMONE LH)

Product code: SC-1460-20




A Fluorescence Immunoassay test kit to detect Luteinizing Hormone (LH) in whole blood, serum or plasma, with the use of Fluorescence Immunoassay Analyzer. The LH Test Cassette (whole blood/serum/plasma) is intended for in vitro quantitative determination Luteinizing Hormone (LH) as an aid in the detection of ovulation. The LH Test Cassette (whole blood/serum/plasma) detects the LH surge signaling that ovulation is likely to occur in the next 24-36 hours. The test utilizes a combination of antibodies including a monoclonal LH antibody to selectively detect elevated levels of LH. Important: The LH surge and ovulation may not occur in all menstrual cycles. For professional in vitro diagnostic use only.

# SCREEN IFA TEST MIOGLOBINA

(FLUORESCENCE IMMUNOASSAY MYOGLOBIN)

Product code: SC-1484-20



A Fluorescence Immunoassay test kit for the diagnosis of Myocardial Infarction (MI) by measuring Myoglobin (MYO) in whole blood, serum or plasma with the use of fluorescence Immunoassay Analyzer. The Myoglobin Test Cassette (whole blood/serum/plasma) is intended for in vitro quantitative determination of Myoglobin in whole blood, serum or plasma as an aid in the diagnosis of Myocardial Infarction (MI). For professional in vitro diagnostic use only.



# SCREEN IFA TEST FSH

(FLUORESCENCE IMMUNOASSAY FOLLICLE STIMULATING HORMONE (FSH))

Product code: SC-1385-20

A Fluorescence Immunoassay for the detection of Follicle Stimulating Hormone (FSH) in whole blood, serum or plasma with the use of Fluorescence Immunoassay Analyzer. The FSH Test Cassette (whole blood/serum/plasma) is intended for in vitro quantitative determination follicle stimulating hormone (FSH) in whole blood, serum or plasma as an aid in the diagnosis of menopause. For professional in vitro diagnostic use only.



# SCREEN IFA TEST INFLUENZA A/B

(FLUORESCENCE IMMUNOASSAY INFLUENZA A/B)

Product code: SC-1408-20

A Fluorescence Immunoassay test kit for the diagnosis of influenza A and influenza B virus in nasopharyngeal swab, throat swab or nasal aspirate specimens with the use of Fluorescence Immunoassay Analyzer. The Influenza A+B Test Cassette (swab/nasal aspirate) is intended for in vitro detection of influenza A and B antigens in nasopharyngeal swab, throat swab or nasal aspirate specimens. It is intended to aid in the rapid differential diagnosis of influenza A and B viral infections. For professional in vitro diagnostic use only.





# SCREEN IFA TEST FOB

(FLUORESCENCE IMMUNOASSAY FECAL OCCULT BLOOD)

Product code: SC-1316-20



A Fluorescence Immunoassay test kit for the quantitative detection of Fecal Occult Blood (FOB) in feces with the use of Fluorescence Immunoassay Analyzer. In the early stages, gastrointestinal problems such as colon cancer, ulcers, polyps, colitis, diverticulitis and fissures may not show any visible symptoms, only occult blood. Traditional guaiac-based methods lack sensitivity and specificity and have diet restrictions prior to testing. The test uses a double antibody sandwich assay to selectively detect Fecal Occult Blood at 10-1000ng/ml. For professional in vitro diagnostic use only.

# SCREEN IFA TEST D-DIMERO

(FLUORESCENCE IMMUNOASSAY D-DIMER)

Product code: SC-1583-20




A Fluorescence Immunoassay test kit for measuring D-Dimer in whole blood or plasma with the use of Fluorescence Immunoassay Analyzer. The D-Dimer Test Cassette (whole blood/plasma) uses Immunofluorescence to measure D-Dimer in whole blood or plasma as an aid in the diagnosis of DVT and PE. The sample moves through the strip from sample pad to absorbent pad. If the specimen contains D-Dimer, it attaches to the fluorescent microspheres-conjugated anti-D-Dimer antibodies. For professional in vitro diagnostic use only.

# SCREEN IFA TEST PROCALCITONINA

(FLUORESCENCE IMMUNOASSAY PROCALCITONIN)

Product code: SC-1392-20

A Fluorescence Immunoassay test kit for the diagnosis of inflammatory condition by measuring Procalcitonin (PCT) in whole blood, serum or plasma with the use of Fluorescence Immunoassay Analyzer. The PCT Test Cassette (whole blood/serum/plasma) is based on Fluorescence immunoassay for the quantitative detection of human Procalcitonin in whole blood, serum or plasma as an aid in the diagnosis of inflammatory conditions. For professional in vitro diagnostic use only.



# SCREEN IFA TEST TSH

(FLUORESCENCE IMMUNOASSAY THYROID STIMULATING HORMONE)

Product code: SC-1378-20

A Fluorescence Immunoassay test kit to detect TSH quantity levels in serum or plasma with the use of Fluorescence Immunoassay Analyzer. TSH Test Cassette (Immunofluorescence Assay) is intended for in vitro quantitative determination of Thyroid Stimulating Hormone (TSH) in serum or plasma. Measurement of TSH is useful to aid in the screening the adult population for primary hypothyroidism by medical professionals. It could also be used in screening neonates for hypothyroidism. For professional in vitro diagnostic use only.





# SCREEN IFA TEST TROPONIN I

(FLUORESCENCE IMMUNOASSAY TROPONIN I)

Product code: SC-1590-20



A Fluorescence Immunoassay test kit for the diagnosis of Myocardial Infarction (MI) to detect Cardiac Troponin I (cTnI) in whole blood, serum or plasma with the use of Fluorescence Immunoassay Analyzer. The Cardiac Troponin I Test Cassette (whole blood/serum/plasma) is intended for in vitro quantitative determination of human Cardiac Troponin I in whole blood, serum or plasma as an aid in the diagnosis of Myocardial Infarction (MI). For professional in vitro diagnostic use only.

# SCREEN IFA TEST CEA

(FLUORESCENCE IMMUNOASSAY CARCINOEMBRYONIC ANTIGEN)

Product code: SC-1354-20



A Fluorescence Immunoassay test kit for measuring Carcinoembryonic Antigen (CEA) in serum or plasma with the use of Fluorescence Immunoassay Analyzer. Carcinoembryonary Antigen (CEA) is an antigen associated with cancer characterized as oncofetal glycoprotein. CEA is expressed in a variety of malignant tumors, in particular lung or gastrointestinal cancers (e.g. colon cancer, liver cancer and lung cancer). For professional in vitro diagnostic use only.



# • PREGNANCY GINECOLOGY •

43

PREGNANCY/GYNECOLOGY

# SCREEN TEST PH VAGINALE

(VAGINAL PH RAPID TEST)

Product code: SC-0753-25



The test is a screening device for use in the detection of the major form of vaginitis in vaginal fluid specimens from women concerned about their vaginal health. Vaginitis is defined as irritation of the vagina, a troublesome condition that affects millions of women of all ages in all parts of the World. For professional in vitro diagnostic use only.

# SCREEN TEST FAMILY SPERM

(SPERM CONCENTRATION RAPID TEST)

Product code: SC-0784-40



The kit is used for qualitatively detecting sperm count of human semen in vitro and conducting auxiliary clinical diagnosis of the infertility and self-evaluation male fertility. For professional in vitro diagnostic use only.

# SCREEN TEST HCG

(PREGNANCY RAPID TEST (10NG))

Product code: SC-0425-40

The HCG Pregnancy Rapid Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of Human Chorionic Gonadotropin in urine to aid in the early detection of pregnancy. The test uses two lines to indicate results. The test utilizes a combination of antibodies including a monoclonal HCG antibody to selectively detect elevated levels of HCG. The HCG Rapid Pregnancy Test is a qualitatively identifying test the presence of HCG in a urine sample with a sensitivity of 10mIU/ml. For professional in vitro diagnostic use only.



# SCREEN TEST HCG

(PREGNANCY RAPID TEST (25NG))

Product code: SC-0432-40

The HCG Pregnancy Rapid Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of Human Chorionic Gonadotropin in urine to aid in the early detection of pregnancy. The test uses two lines to indicate results. The test utilizes a combination of antibodies including a monoclonal HCG antibody to selectively detect elevated levels of HCG. The HCG Rapid Pregnancy Test is a qualitatively identifying test the presence of HCG in a urine sample with a sensitivity of 25mIU/ml. For professional in vitro diagnostic use only.





# • INFECTIVE • • DISEASES •

46

INFECTIVE DISEASES

# SCREEN TEST COVID-19 SALIVA

(SARS-COV-2 ANTIGEN ORAL FLUID RAPID TEST)

Product code: SC-1286-20

The Covid-19 Antigen Rapid Test (Oral Fluid) is a rapid chromatographic immunoassay for the qualitative detection of SARS-CoV-2 Antigens in oral fluid specimens from individuals with suspected SARS-CoV-2 infection in conjunction with clinical presentation and the results of other laboratory tests. An antigen is generally detectable in upper respiratory specimens during the acute phase of infection. For professional in vitro diagnostic use only.



# SCREEN TEST COVID-19 RBD-IGG

(SARS-COV-2 RBD-IGG RAPID TEST)

Product code: SC-1415-20

The Covid-19 IgG Rapid Test (whole blood/serum/plasma) is a lateral flow chromatographic immunoassay for the qualitative detection of IgG antibodies to SARS-CoV2 Spike (S) protein Receptor Binding Domain (RBD) and/or Nucleocapsid protein in human whole blood, serum or plasma specimen. The COVID-19 IgG Rapid Test (whole blood/serum/plasma) is intended for use as an aid in identifying individuals with an adaptive immune response to SARS-CoV-2. For professional in vitro diagnostic use only.



# SCREEN TEST COVID-19 + FLU A/B

(SARS-COV-2 ANTIGEN + FLU A/B RAPID TEST)

Product code: SC-1262-20



A rapid chromatographic immunoassay for the qualitative detection of Sars-CoV-2, Influenza A and Influenza B virus antigens present in human nasopharynx. An antigen is generally detectable in upper respiratory specimens during the acute phase of infection. Positive results indicate the presence of viral antigens but clinical correlation with patient history and other diagnostic information are necessary to determine infection status. For professional in vitro diagnostic use only.

# SCREEN TEST COVID-19 IGG/IGM

(SARS-COV-2 IGG/IGM RAPID TEST)

Product code: SC-1200-10



A rapid test for the qualitative detection of IgG and IgM antibodies to 2019-nCov, Sars-CoV-2 in human Fingerstick Whole Blood specimens. The Sars-CoV-2, 2019-nCOV IgG/IgM Rapid Test Cassette is a lateral flow chromatographic immunoassay for the qualitative detection of IgG and IgM antibodies to 2019-nCov in human whole blood, serum or plasma specimen. For professional in vitro diagnostic use only.



# SCREEN TEST COVID-19 AG

(SARS-COV-2 ANTIGEN RAPID TEST)

Product code: SC-1255-20

A rapid chromatographic immunoassay for the qualitative detection of Sars-CoV-2 antigens present in human nasopharynx. An antigen is generally detectable in upper respiratory specimens during the acute phase of infection. Positive results indicate the presence of viral antigens but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease. For professional in vitro diagnostic use only.



# SCREEN TEST COVID-19 NT

(SARS-COV-2 NEUTRALIZATION ANTIBODY RAPID TEST)

Product code: SC-1569-20

A rapid chromatographic immunoassay for the qualitative detection of Neutralization Antibodies against Sars-CoV-2 that block the interaction between the receptor binding domain of the viral spike Glycoprotein (RBD) with the cell surface receptor ACE2 in human serum, plasma or whole blood. It is intended for use as an aid in identifying individuals with an adaptive immune response to Sars-CoV-2. This test should not be used to diagnose acute Sars-CoV-2 infection. For professional in vitro diagnostic use only.



# SCREEN TEST GIARDIA LAMBLIA

(GIARDIA LAMBLIA RAPID TEST)

Product code: SC-0555-20



A rapid chromatographic immunoassay for the qualitative detection of Giardia Lamblia antigen in human feces specimen. Parasitic infections remain a very serious health problem worldwide. Giardia Lamblia is the most common protozoa known to be responsible for one of the main causes of severe diarrhea in humans particularly in immunocompromised people. For professional in vitro diagnostic use only.

# SCREEN TEST ROTAVIRUS

(ROTAVIRUS RAPID TEST)

Product code: SC-0692-25



A rapid one-step test for the qualitative detection of Rotavirus in human feces. The Rotavirus Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of Rotavirus in human feces specimens to aid in the diagnosis of Rotavirus infection. For professional in vitro diagnostic use only.

# SCREEN TEST CLOSTRIDIUM A+B

(CLOSTRIDIUM DIFFICILE TOXIN A+TOXIN B RAPID TEST)

Product code: SC-0494-10

A rapid chromatographic immunoassay for the qualitative detection of Clostridium Difficile Toxin A and Toxin B antigens in the human feces specimen. Clostridium Difficile Toxin A + Toxin B Combo Rapid Test Cassette detects two distinct antigens in fecal specimens for C. Difficile, Viz., Toxin A and Toxin B on two different test strips in a single test cassette, thus simultaneously detecting two antigens specific to Clostridium difficile. For professional in vitro diagnostic use only.



# SCREEN TEST H.PYLORI AB

(H. PYLORI ANTIBODY RAPID TEST)

Product code: SC-0609-40

A rapid test for the qualitative detection of Antibodies to Helicobacter Pylori (H. Pylori) in whole blood, serum or plasma. The H. Pylori Antibody Rapid Test Cassette (whole blood/serum/plasma) is a rapid chromatographic immunoassay for the qualitative detection of Antibodies to H. Pylori in whole blood, serum or plasma to aid in the diagnosis of H. Pylori Infection. For professional in vitro diagnostic use only.





# SCREEN TEST CHLAMYDIA

(CHLAMYDIA RAPID TEST)

Product code: SC-0500-20



A rapid test for the qualitative detection of Chlamydia antigen in female cervical swab, male urethral swab and male urine specimens. The Chlamydia Rapid Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of Chlamydia trachomatis in female cervical swab, male urethral swab and male urine specimens to aid in the diagnosis of Chlamydia infection. For professional in vitro diagnostic use only.

# SCREEN TEST CLOSTRIDIUM

(CLOSTRIDIUM DIFFICILE GDH RAPID TEST)

Product code: SC-0470-25



A rapid chromatographic immunoassay for the qualitative detection of Clostridium Difficile GDH in the human feces specimen. This is a ready-to-use test that is based on the use of a membrane technology with colloidal gold. A nitrocellulose membrane is sensitized with antibody directed against Clostridium Difficile Antigen (GDH). For professional in vitro diagnostic use only.

# SCREEN TEST MONONUCLEOSIS

(MONONUCLEOSIS RAPID TEST)

Product code: SC-0647-25

A rapid test for the diagnosis of Infectious Mononucleosis (IM) to detect Infectious Mononucleosis heterophile antibodies qualitatively in whole blood, serum or plasma. The MONO Rapid Test Cassette (whole blood/serum/plasma) is a simple test that utilizes an extract of bovine erythrocytes to qualitatively and selectively detect Infectious Mononucleosis heterophile antibodies in whole blood, serum or plasma in minutes. For professional in vitro diagnostic use only.



# SCREEN TEST MYCOPLASMA

(MYCOPLASMA PNEUMONIAE ANTIGEN TEST)

Product code: SC-0654-25

Mycoplasma Pneumoniae Antigen Rapid Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of Mycoplasma Pneumoniae (M. Pneumoniae) antigens in human throat swabs. It is intended to aid in the rapid differential diagnosis of Mycoplasma Pneumoniae infections. For professional in vitro diagnostic use only.



# SCREEN TEST H.PYLORI AG

(H.PYLORI ANTIGEN RAPID TEST)

Product code: SC-0616-25




A rapid test for the qualitative detection of Helicobacter Pylori (H. Pylori) Antigens in human feces. The H. Pylori Antigen Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of H. Pylori Antigens in human feces specimens to aid in the diagnosis of H. Pylori infection. For professional in vitro diagnostic use only.

# SCREEN TEST FLU A/B

(INFLUENZA A+B RAPID TEST)

Product code: SC-0623-20



A rapid chromatographic immunoassay for the qualitative detection of influenza virus, including the Swine flu H1N1 H3N2. The result is divided into Influenza A (including H1N1 and H3N2) and B type (less dangerous). Available for: hospitals, government bodies, physicians, analysis labs, companies, retailers, etc. For professional in vitro diagnostic use only.



(STREPTOCOCCUS A RAPID TEST)

Single-phase immunoassay on a membrane for the rapid determination and qualitative detection of Antigens of Group A Streptococcus directly from the pharyngeal swab also included in the package the positive and negative control. For professional in vitro diagnostic use only.



(TETANUS RAPID TEST)

Rapid immunoassay for the qualitative detection of Tetanus Antitoxin Antibodies on blood samples. The test is designed to assist tetanus prophylaxis of patients at risk. The immunoprotection is given by the measurement of the concentration of anti-tetanus antibodies. As part of emergency medicine, to know in real time the immune status of patients at risk (deep wounds) allows a better management of the appropriate Tetanus prophylaxis. For professional in vitro diagnostic use only.



# SCREEN TEST NOROVIRUS

(NOROVIRUS RAPID TEST)

Product code: SC-0678-25



The Norovirus Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of Norovirus in human feces specimens to aid in the diagnosis of Norovirus infection. For professional in vitro diagnostic use only.

# SCREEN TEST SIFILIDE

(SYPHILIS RAPID TEST)

Product code: SC-0722-40



The Syphilis Rapid Test Cassette (whole blood/serum/plasma) is a rapid chromatographic immunoassay for the qualitative detection of antibodies (IgG and IgM) to Treponema Pallidum (TP) in whole blood, serum or plasma to aid in the diagnosis of Syphilis. The Syphilis Rapid Test Cassette (whole blood/serum/plasma) is a qualitative membrane-based immunoassay for the detection of TP antibodies (IgG and IgM) in whole blood, serum or plasma. For professional in vitro diagnostic use only.

# SCREEN TEST RUBOLIA IGM

(RUBELLA IGM RAPID TEST)

Product code: SC-0685-40

A rapid test for the qualitative detection of Rubella IgM antibodies in serum or plasma. The Rubella IgM Rapid Test Cassette is a lateral flow chromatographic immunoassay for the qualitative detection of Rubella IgM antibodies in serum or plasma to aid in the diagnosis of Rubella infection. For professional in vitro diagnostic use only.



# SCREEN TEST LEGIONELLA

(LEGIONELLA PNEUMOPHILA RAPID TEST)

Product code: SC-0630-25

Legionella Pneumophila Rapid Test Cassette (Urine) is an in vitro diagnostic test based on immunochromatographic assay. It is designed for detection of soluble antigen from Legionella Pneumophila serogroup 1 in human urine specimen. For professional in vitro diagnostic use only.





# SCREEN TEST CHIKUNGUNYA

(CHIKUNGUNYA IGG/IGM RAPID TEST)

Product code: SC-0517-20



The Chikungunya IgG/IgM Rapid Test Cassette (whole blood/serum/plasma) is a rapid chromatographic immunoassay for the qualitative detection of IgG and IgM antibodies to Chikungunya in human's whole blood, serum or plasma. It is intended to be used as a screening test and as an aid in the diagnosis of infection with CHIKV. For professional in vitro diagnostic use only.

# SCREEN TEST ZIKA

(NS1 ZIKA RAPID TEST)

Product code: SC-0760-40



The NS1 Zika Rapid Test Cassette (whole blood/serum/plasma) is a rapid chromatographic immunoassay for the qualitative detection of NS1 antibodies to Zika in human's whole blood, serum or plasma. It is intended to be used as a screening test and as an aid in the diagnosis of infection with Zika. For professional in vitro diagnostic use only.

## (MALARIA P.F. /P.V. RAPID TEST)

The Malaria P.f./P.v. Rapid Test Cassette (whole blood) is a rapid chromatographic immunoassay for the qualitative detection of two kinds of circulating Plasmodium Falciparum (P. falciparum (P.f.) and P. vivax (P.v.)) in whole blood. For professional in vitro diagnostic use only.



(ADENOVIRUS RAPID TEST)

A rapid one step test for the qualitative detection of Adenovirus in human feces. Rapid and accurate diagnosis of gastroenteritis due to adenovirus is helpful in establishing the etiology of gastroenteritis and related patient management. The Adenovirus Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of Adenovirus in human feces specimens to aid in the diagnosis of adenovirus infection. For professional in vitro diagnostic use only.



# SCREEN TEST PNEUMOCOCCO

(S. PNEUMONIAE ANTIGEN RAPID TEST)

Product code: SC-0708-25



The S. Pneumoniae Antigen Rapid Test Cassette (urine) is a rapid chromatographic immunoassay for the qualitative detection of Streptococcus Pneumoniae Antigens in human urine specimen. For professional in vitro diagnostic use only.

# SCREEN TEST TUBERCULOSI

(TUBERCULOSIS RAPID TEST)

Codice Prodotto: SC-0739-40



Rapid test for the qualitative detection of Anti-TB antibodies (IgG, IgM and IgA isotypes) in whole blood, serum or plasma samples. The Rapid Tuberculosis Test (whole blood/serum/plasma) is a qualitative, solid phase, direct, two-site immunoassay for the detection of Anti-TB antibodies in whole blood, serum or plasma samples. For professional in vitro diagnostic medical use only.



# SCREEN TEST DENGUE NS1

(DENGUE NS1 RAPID TEST)

Product code: SC-0548-25

A rapid chromatographic immunoassay for the qualitative detection of NS1 antigen of Dengue virus in human whole blood, serum or plasma as an aid in the diagnosis of Dengue infections. For professional in vitro diagnostic use only.



# SCREEN TEST DENGUE COMBO

(DENGUE COMBO RAPID TEST)

Product code: SC-0531-25

A rapid chromatographic immunoassay for the qualitative detection of NS1 antigen and IgG and IgM antibodies of Dengue virus in human whole blood, serum or plasma as an aid in the diagnosis of Dengue infections. For professional in vitro diagnostic use only.



# SCREEN TEST GONORRHEA

(GONORRHEA RAPID TEST)

Product code: SC-1187-25



A rapid test for the qualitative detection of Gonorrhea antigen in female cervical swab and male urethral swab specimens. The Gonorrhea Rapid Test Cassette (swab) is a rapid chromatographic immunoassay for the qualitative detection of *Neisseria Gonorrhea* in female cervical swab and male urethral swab specimens to aid in the diagnosis of Gonorrhea infection. Gonorrhea is a sexually transmitted disease caused by the bacterium *Neisseria Gonorrhoeae*. For professional in vitro diagnostic use only.

# SCREEN TEST CRYPTOSPORIDIUM

(CRYPTOSPORIDIUM ANTIGEN TEST )

Product code SC-0524-10



A rapid chromatographic immunoassay for the qualitative detection of *Cryptosporidium* Antigens in human feces. Cryptosporidiosis is a diarrhoeal disease caused by microscopic parasites of the genus *Cryptosporidium*. Once an animal or person is infected, the parasite lives in the intestine and passes in the stool. For professional in vitro diagnostic use only.

# SCREEN TEST EPATITE A - HAV

(HAV IGM RAPID TEST)

Product code: SC-0562-25

The HAV IgM Rapid Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of IgM antibody to Hepatitis A Virus (HAV) in serum or plasma. For professional in vitro diagnostic use only.



# SCREEN TEST EPATITE B - HBsAg

(HBsAg RAPID TEST)

Product code: SC-0579-20

The HBsAg Rapid Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of Hepatitis B Surface Antigen in whole blood, serum or plasma. Viral hepatitis is a systemic disease primarily involving the liver. Most cases of acute viral hepatitis are caused by Hepatitis A Virus, Hepatitis B Virus (HBV) or Hepatitis C Virus. The complex antigen found on the surface of HBV is called HBsAg. Previous designations included the Australia or Au antigen. For professional in vitro diagnostic use only.





# SCREEN TEST EPATITE C - HCV

(HCV RAPID TEST)

Product code: SC-0586-20



A rapid chromatographic immunoassay for the qualitative detection of antibody to Hepatitis C Virus in whole blood, serum or plasma. Hepatitis C Virus (HCV) is a small, enveloped, positive-sense, single-stranded RNA Virus. HCV is now known to be the major cause of parenterally transmitted non-A, non-B hepatitis. For professional in vitro diagnostic use only.

# SCREEN TEST CAMPYLOBACTER

(CAMPYLOBACTER RAPID TEST)

Product code: SC-0463-10



A rapid chromatographic immunoassay for the qualitative detection of Campylobacter antigen in human feces specimen. The Campylobacter Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of Campylobacter antigen in human feces specimen. For professional in vitro diagnostic use only.

# SCREEN TEST HIV 1.2

(HIV 1.2 RAPID TEST)

Product code: SC-0593-20

HIV Rapid Test (whole blood/serum/plasma) is an immunochromatographic test for the detection qualitative antibody to Human Immunodeficiency Virus (HIV) in whole blood, serum or plasma to help diagnose HIV infection. The HIV Rapid Test Cassette (whole blood/serum/plasma) is a qualitative immunoassay, membrane-based for the detection of HIV antibodies in whole blood, serum or plasma. For professional in vitro diagnostic use only.



# INFLAMMATORY ... MARKERS ...

66

INFLAMMATORY MARKERS



# SCREEN TEST CRP

(C-REACTIVE PROTEIN RAPID TEST)

Product code: SC-0975-25

A rapid chromatographic immunoassay for the semi-quantitative result of CRP levels in whole blood, serum or plasma. The CRP Rapid Test Cassette (whole blood/serum/plasma) is a rapid chromatographic immunoassay for the semi-quantitative detection of human CRP in whole blood, serum or plasma as an aid in the diagnosis of Myocardial Infarction (MI). For professional in vitro diagnostic use only



# SCREEN TEST PROCALCITONINA

(PROCALCITONIN (PCT) RAPID TEST)

Product code: SC-0999-25

The PCT Rapid Test Cassette (Whole Blood/Serum/Plasma) is a rapid chromatographic immunoassay for the qualitative detection of Procalcitonin in serum or plasma. The diagnostic value of PCT is important due to the close correlation between PCT concentration and the severity of inflammation. For professional in vitro diagnostic use only.



# ...TUMORAL MARKERS...

68

TUMORAL MARKERS

# SCREEN TEST FOB

(FECAL OCCULT BLOOD RAPID TEST)

Product code: SC-0883-25

A rapid one step test for the qualitative detection of Human Occult Blood in feces. The FOB Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of human occult blood in feces. Many diseases can cause hidden blood in the feces. This is also known as Fecal Occult Blood (FOB), Human Occult Blood or Human Hemoglobin. In the early stages, gastrointestinal problems such as colon cancer, ulcers, polyps, colitis, diverticulitis and fissures may not show any visible symptoms only occult blood. For professional in vitro diagnostic use only.



# SCREEN TEST PSA

(PROSTATE-SPECIFIC ANTIGEN RAPID TEST)

Product Code: SC-0890-40

A rapid test for the qualitative detection of Prostate Specific Antigen (PSA) in whole blood, serum or plasma. The PSA Prostate Specific Antigen Semi-Quantitative Rapid Test Cassette (whole blood /serum/plasma) is a rapid chromatographic immunoassay for semi-quantitative detection of Prostate Specific Antigen in whole blood, serum or plasma. For professional in vitro diagnostic use only.





# SCREEN TEST CALPROTECTINA

(CALPROTECTIN RAPID TEST)

Product Code: SC-0777-25



The Calprotectin Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of Calprotectin in human feces specimen which might be useful for the diagnosis of inflammatory gastrointestinal disorders. The Calprotectin Rapid Test Cassette (Feces) can detect levels of Calprotectin as low as 50  $\mu\text{g/g}$  or 140  $\text{ng/ml}$  feces. For professional in vitro diagnostic use only.

# ...CARDIAC... MARKERS.

71

CARDIAC MARKERS

# SCREEN TEST MIOG./TROP./CK-MB

(MYOGLOBIN/CK-MB/TROPONIN I COMBO TEST)

Product code: SC-0982-40



A rapid chromatographic immunoassay to detect Myoglobin, CK-MB and cardiac Troponin I (cTnI) qualitatively in whole blood, serum or plasma. The Myoglobin/CK-MB/Troponin I Combo Rapid Test Cassette (whole blood/serum/plasma) is a rapid chromatographic immunoassay of support in the diagnosis of Myocardial Infarction (MI). For professional in vitro diagnostic use only.



# ..URINE.. ANALYSIS

**73** URINE ANALYSIS

# SCREEN GENERIC TEST URINA 10

(URINANALYSIS REAGENT STRIPS)

Product code: SC-0906-100



Reagent Strips for rapid detection of multiple analytes in human urine. The Urinalysis Reagent Strips (Urine) are firm plastic strips onto which several separate reagent areas are affixed. The test is for the qualitative and semi-quantitative detection of one or more of the following analytes in urine: Leukocytes, Nitrite, Urobilinogen, Protein, PH, Blood, Specific Gravity, Ketone, Bilirubin, Glucose. The Urinalysis Reagent Strips (Urine) are for single use in professional near-patient (point-of-care) and centralized laboratory locations. For professional in vitro diagnostic use only.

# SCREEN GENERIC TEST URINA 11

(URINANALYSIS REAGENT STRIPS)

Product code: SC-0913-100



Reagent Strips for rapid detection of multiple analytes in human urine. The Urinalysis Reagent Strips (Urine) are firm plastic strips onto which several separate reagent areas are affixed. The test is for the qualitative and semi-quantitative detection of one or more of the following analytes in urine: Ascorbic acid, Glucose, Bilirubin, Ketone, Specific Gravity, Blood, PH, Protein, Urobilinogen, Nitrite, Leukocytes. The Urinalysis Reagent Strips (Urine) are for single use in professional near-patient (point-of-care) and centralized laboratory locations. For professional in vitro diagnostic use only.

# SCREEN GENERIC TEST URINA 12

(URINANALYSIS REAGENT STRIPS)

Product code: SC-0920-100

Reagent Strips for rapid detection of multiple analytes in human urine. The Urinalysis Reagent Strips (Urine) are firm plastic strips onto which several separate reagent areas are affixed. The test is for the qualitative and semi-quantitative detection of one or more of the following analytes in urine: Ascorbic Acid, Glucose, Bilirubin, Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrites, Leukocytes and Microalbuminuria. The Urinalysis Reagent Strips (Urine) are for single use in professional near-patient (point-of-care) and centralized laboratory locations. For professional in vitro diagnostic use only.





# ...OTHER TESTS...

# SCREEN TEST VITAMINA D

(VITAMIN D RAPID TEST)

Product code: SC-0791-25

A rapid chromatographic immunoassay for the semi-quantitative detection of 25-Hydroxyvitamin D (25 (OH) D) in human fingerstick whole blood at a cut-off concentration of  $30 \pm 4\text{ng/mL}$ . Vitamin D refers to a group of fat-soluble secosteroids responsible for increasing intestinal absorption of calcium, iron, magnesium, phosphate, and zinc. Vitamin D is transported to the liver where it is metabolized to 25-Hydroxy Vitamin D. The Vitamin D Rapid Test Cassette provides a preliminary diagnostic test result and can be used to screening for Vitamin D deficiency. For professional in vitro diagnostic use only.



# RAPID TEST READERS

78

RAPID TEST READERS



# SCREEN ANALYSIS READER

(URINE ANALYSIS READER)

Product code: SC-0937-01 10 parameters / SC-0944-01 11 parameters / SC-0951-01 12 parameters

Reader for Generic Urine Reagent Strips of 10, 11 and 12 Parameters with Thermal Printer. Small size analyzer, accurate and effective, maintenance-free, wireless thermal printer (optional), capacitive touch keyboard, sensitive reaction, internal rechargeable (1200mAh battery).



# SCREEN READER PLUS

(RAPID TEST READER TOUCH SCREEN)

Product Code: SC-1750-01

Instrument composed of a PC with integrated Screen Rapid Test Reader. A single tool that reads and interprets Urine Drug Test Results Immersion / Cassette, Saliva Drug Test, Tetanus Test, Tuberculosis Test, Syphilis Test, Streptococcus Test, PCT Test, Mycoplasma Test, H. Pylori Test, Fob Test, PCR Test, Calprotectin Test, Alcohol Test Etc, Hcg Test, Covid 19 Test, etc. The tool allows you to carry out the objective print of the test result with the possibility of archiving the results in Pdf files for your own archive.



# SCREEN READER SOFTWARE + SCREEN BOX

(SOFTWARE) (RAPID TEST READER)

Product code: SC-1194-01 / SC-0869-01



Screen Reader is a unique software that allows reading and reporting objectively the rapid test results, which can be printed or archived as PDF files. Screen Reader Software is compatible with all windows operating systems. A unique tool that reads and interprets the results of Urine or Saliva drug tests, Tetanus test, Streptococcus test, Tuberculosis test, Fob test, single-use ET G urine alcohol test, Norovirus test, Mycoplasma Test, HCG Pregnancy Test, PCR Test, Covid-19 rapid test, etc.





**SCREEN ITALIA S.r.l.**

Via dell'Artigianato, 16  
06089 – Torgiano  
Perugia – Italia  
VAT: IT03137260547  
Tel. +39 388 1780178  
Whatsapp +39 388 1780178  
Skype sales@screenitalia.it  
sales@screenitalia.it  
www.screenitalia.eu